

# Dispensette® S Bottletop Dispenser

Setting the standard  
for over a half century



**BRANDTECH**<sup>®</sup>  
SCIENTIFIC, INC.

# BRAND® Dispensette® S Bottletop Dispensers

## Easy calibration technique

Calibration and adjustments according to ISO 9001 and GLP are done within seconds

## Calibration mechanism

Visible flag indicates calibration has been adjusted from factory specifications

## Discharge valve with safety ball

Closes when discharge tube is not mounted to prevent accidental dispensing

## Large viewing window

Allows easy priming verification

## Volume selection

With interior scalloped track for reproducible volume setting

## Recirculation Valve

Eliminates reagent waste and splashing during priming

## Freely rotating valve block

GL 45 mm thread allows bottle label to always face user

## Olive-shaped filling valve

For firm filling tube attachment

## Hinged screw cap

Swings out of the way when dispensing

## Threaded safety cap


Coarse thread allows fingertip on/off

## Recirculation tube

Safely returns reagent to the bottle

## Telescoping filling tube

Adjusts easily to a broad range of bottle sizes. No measuring or cutting is required



Dispensette® S Organic Digital  
without recirculation valve

Dispensette® S Analog-adjustable with  
recirculation valve

# BRAND® Dispensette® S Bottletop Dispensers

Dispensette® S						Without recirculation valve		With recirculation valve		
Volume, mL	Increments, mL	A* < ±		CV* ≤		Cat. No.	2021	Cat. No.	2021	
		%	µL	%	µL		List Price		List Price	
<b>Dispensette® S, Digital</b>										
0.1-1	0.005	0.6	6	0.2	2	<a href="#">4600310</a>	\$534.80	<a href="#">4600311</a>	\$568.00	
0.2-2	0.01	0.5	10	0.1	2	<a href="#">4600320</a>	534.80	<a href="#">4600321</a>	568.00	
0.5-5	0.02	0.5	25	0.1	5	<a href="#">4600330</a>	534.80	<a href="#">4600331</a>	568.00	
1-10	0.05	0.5	50	0.1	10	<a href="#">4600340</a>	534.80	<a href="#">4600341</a>	568.00	
2.5-25	0.1	0.5	125	0.1	25	<a href="#">4600350</a>	722.20	<a href="#">4600351</a>	755.40	
5-50	0.2	0.5	250	0.1	50	<a href="#">4600360</a>	732.80	<a href="#">4600361</a>	766.00	
<b>Dispensette® S, Analog-adjustable</b>										
0.1-1	0.02	0.6	6	0.2	2	<a href="#">4600100</a>	468.80	<a href="#">4600101</a>	490.00	
0.2-2	0.05	0.5	10	0.1	2	<a href="#">4600120</a>	468.80	<a href="#">4600121</a>	490.00	
0.5-5	0.1	0.5	25	0.1	5	<a href="#">4600130</a>	468.80	<a href="#">4600131</a>	490.00	
1-10	0.2	0.5	50	0.1	10	<a href="#">4600140</a>	468.80	<a href="#">4600141</a>	490.00	
2.5-25	0.5	0.5	125	0.1	25	<a href="#">4600150</a>	656.20	<a href="#">4600151</a>	677.60	
5-50	1.0	0.5	250	0.1	50	<a href="#">4600160</a>	672.20	<a href="#">4600161</a>	694.40	
10-100	1.0	0.5	500	0.1	100	<a href="#">4600170</a>	1,047.20	<a href="#">4600171</a>	1,069.60	
<b>Dispensette® S, Fixed-volume</b>										
1		0.6	6	0.2	2	<a href="#">4600210</a>	468.80	<a href="#">4600211</a>	490.00	
2		0.5	10	0.1	2	<a href="#">4600220</a>	468.80	<a href="#">4600221</a>	490.00	
5		0.5	25	0.1	5	<a href="#">4600230</a>	468.80	<a href="#">4600231</a>	490.00	
10		0.5	50	0.1	10	<a href="#">4600240</a>	468.80	<a href="#">4600241</a>	490.00	
<b>Dispensette® S Organic</b>										
Dispensette® S Organic						Without recirculation valve		With recirculation valve		
Volume, mL	Increments, mL	A* < ±		CV* ≤		Cat. No.	2021	Cat. No.	2021	
		%	µL	%	µL		List Price		List Price	
<b>Dispensette® S Organic, Digital</b>										
0.5-5	0.02	0.5	25	0.1	5	<a href="#">4630330</a>	\$589.40	<a href="#">4630331</a>	\$622.20	
1-10	0.05	0.5	50	0.1	10	<a href="#">4630340</a>	589.40	<a href="#">4630341</a>	622.20	
2.5-25	0.1	0.5	125	0.1	25	<a href="#">4630350</a>	793.80	<a href="#">4630351</a>	826.60	
5-50	0.2	0.5	250	0.1	50	<a href="#">4630360</a>	809.60	<a href="#">4630361</a>	842.60	
<b>Dispensette® S Organic, Analog-adjustable</b>										
0.5-5	0.1	0.5	25	0.1	5	<a href="#">4630130</a>	507.00	<a href="#">4630131</a>	540.20	
1-10	0.2	0.5	50	0.1	10	<a href="#">4630140</a>	507.00	<a href="#">4630141</a>	540.20	
2.5-25	0.5	0.5	125	0.1	25	<a href="#">4630150</a>	710.40	<a href="#">4630151</a>	743.60	
5-50	1.0	0.5	250	0.1	50	<a href="#">4630160</a>	727.80	<a href="#">4630161</a>	760.60	
10-100	1.0	0.5	500	0.1	100	<a href="#">4630170</a>	1,141.00	<a href="#">4630171</a>	1,174.20	
<b>Dispensette® S Organic, Fixed-volume</b>										
5		0.5	25	0.1	5	<a href="#">4630230</a>	507.00	<a href="#">4630231</a>	540.20	
10		0.5	50	0.1	10	<a href="#">4630240</a>	507.00	<a href="#">4630241</a>	540.20	



Dispensette® S



Dispensette® S Organic



Dispensette® S Trace Analysis

Dispensette® S Trace Analysis						Without recirculation valve		With recirculation valve		
Volume, mL	Valve Spring	A* < ±		CV* ≤		Cat. No.	2020	Cat. No.	2020	
		%	µL	%	µL		List Price		List Price	
<b>Dispensette® S Trace Analysis, Analog-adjustable</b>										
1-10	Platinum-iridium	0.5	50	0.1	10	<a href="#">4640040</a>	\$1,151.40	<a href="#">4640041</a>	\$1,183.80	
1-10	Tantalum	0.5	50	0.1	10	<a href="#">4640240</a>	1,151.40	<a href="#">4640241</a>	1,183.80	

A\*=Accuracy, CV\*=Coefficient of Variation

\* The value of accuracy and coefficient of variation are final test values referring to the delivered nominal volume, instrument and distilled water at equilibrium with ambient temperature (20 °C/68 °F) and with smooth operation.



**Flexible Discharge Tube**  
(For more information visit [www.brandtech.com](http://www.brandtech.com))

# Areas of Application / Suggested Dispenser

## ■ Dispensette® S (Disp. S) ■ Dispensette® S Organic (Disp. S Organic)

Reagent	Disp. S	Disp. S Organic	Reagent	Disp. S	Disp. S Organic	Reagent	Disp. S	Disp. S Organic
Acetaldehyde	+	+	Cyclohexane		+	Methylene chloride		+
Acetic acid (glacial), 100%	+	+	Cyclohexanone	+	+	Mineral oil (Engine oil)	+	+
Acetic acid, ≤ 96%	+	+	Cyclopentane		+	Monochloroacetic acid	+	+
Acetic anhydride		+	Decane	+	+	Nitric acid, ≤ 30%	+	+
Acetone	+	+	1-Decanol	+	+	Nitric acid, 30-70% */ **		+
Acetonitrile	+	+	Dibenzyl ether	+	+	Nitrobenzene	+	+
Acetophenone		+	Dichloroacetic acid		+	Oleic acid	+	+
Acetyl chloride		+	Dichlorobenzene	+	+	Oxalic acid	+	
Acetylacetone	+	+	Dichloroethane		+	n-Pentane		+
Acrylic acid	+	+	Dichloroethylene		+	Peracetic acid		+
Acrylonitrile	+	+	Dichloromethane		+	Perchloric acid	+	+
Adipic acid	+		Diesel oil (Heating oil), bp 250-350° C		+	Perchloroethylene		+
Allyl alcohol	+	+	Diethanolamine	+	+	Petroleum, bp 180-220° C		+
Aluminium chloride	+		Diethyl ether		+	Petroleum ether, bp 40-70° C		+
Amino acids	+		Diethylamine	+	+	Phenol	+	+
Ammonia, ≤ 20%	+	+	1,2 Diethylbenzene	+	+	Phenylethanol	+	+
Ammonia, 20-30%		+	Diethylene glycol	+	+	Phenylhydrazine	+	+
Ammonium chloride	+		Dimethyl sulfoxide (DMSO)	+	+	Phosphoric acid, ≤ 85%	+	+
Ammonium fluoride	+		Dimethylaniline	+	+	Phosphoric acid, 85% + Sulfuric acid, 98%, 1:1	+	+
Ammonium sulfate	+		Dimethylformamide (DMF)	+	+	Piperidine	+	+
n-Amyl acetate	+	+	1,4 Dioxane		+	Potassium chloride	+	
Amyl alcohol (Pentanol)	+	+	Diphenyl ether	+	+	Potassium dichromate	+	
Amyl chloride (Chloropentane)		+	Essential oil		+	Potassium hydroxide	+	
Aniline	+	+	Ethanol	+	+	Potassium permanganate	+	
Barium chloride	+		Ethanolamine	+	+	Propionic acid	+	+
Benzaldehyde	+	+	Ethyl acetate	+	+	Propylene glycol (Propanediol)	+	+
Benzene (Benzol)	+	+	Ethylbenzene		+	Pyridine	+	+
Benzine (Petroleum benzin), bp 70-180° C		+	Ethylene chloride		+	Pyruvic acid	+	+
Benzoyl chloride	+	+	Fluoroacetic acid		+	Salicylaldehyde	+	+
Benzyl alcohol	+	+	Formaldehyde, ≤ 40%	+		Scintillation fluid	+	+
Benzylamine	+	+	Formamide	+	+	Silver acetate	+	
Benzylchloride	+	+	Formic acid, ≤ 100%		+	Silver nitrate	+	
Boric acid, ≤ 10%	+	+	Glycerol	+	+	Sodium acetate	+	
Bromobenzene	+	+	Glycol (Ethylene glycol)	+	+	Sodium chloride	+	
Bromonaphthalene	+	+	Glycolic acid, ≤ 50%	+		Sodium dichromate	+	
Butanediol	+	+	Heating oil (Diesel oil), bp 250-350° C		+	Sodium fluoride	+	
1-Butanol	+	+	Heptane		+	Sodium hydroxide, ≤ 30%	+	
n-Butyl acetate	+	+	Hexane		+	Sodium hypochlorite	+	
Butyl methyl ether	+	+	Hexanoic acid	+	+	Sulfuric acid, ≤ 98%	+	+
Butylamine	+	+	Hexanol	+	+	Tartaric acid	+	
Butyric acid	+	+	Hydriodic acid, ≤ 57% **	+	+	Tetrachloroethylene		+
Calcium carbonate	+		Hydrobromic acid	+	+	Tetrahydrofuran (THF) */ **		+
Calcium chloride	+		Hydrochloric acid, ≤ 20%	+	+	Tetramethylammonium hydroxide	+	
Calcium hydroxide	+		Hydrochloric acid, 20-37% **		+	Toluene		+
Calcium hypochlorite	+		Hydrogen peroxide, ≤ 35%		+	Trichloroacetic acid		+
Carbon tetrachloride		+	Isoamyl alcohol	+	+	Trichlorobenzene		+
Chloro naphthalene	+	+	Isobutanol	+	+	Trichloroethane		+
Chloroacetaldehyde, ≤ 45%	+	+	Isooctane		+	Trichloroethylene		+
Chloroacetic acid	+	+	Isopropanol (2-Propanol)	+	+	Trichlorotrifluoro ethane		+
Chloroacetone	+	+	Isopropyl ether	+	+	Triethanolamine	+	+
Chlorobenzene	+	+	Lactic acid	+	+	Triethylene glycol	+	+
Chlorobutane	+	+	Methanol	+	+	Trifluoro ethane		+
Chloroform		+	Methoxybenzene	+	+	Trifluoroacetic acid (TFA)		+
Chlorosulfonic acid		+	Methyl benzoate	+	+	Turpentine		+
Chromic acid, ≤ 50%	+	+	Methyl butyl ether	+	+	Urea	+	
Chromosulfuric acid	+		Methyl ethyl ketone	+	+	Xylene		+
Copper sulfate	+		Methyl formate	+	+	Zinc chloride, ≤ 10%	+	
Cresol		+	Methyl propyl ketone	+	+	Zinc sulfate, ≤ 10%	+	
Cumene (Isopropyl benzene)	+	+						

\* Choose ETFE/PTFE adapters, if required, \*\* use PTFE seal for valve block catalog number #704486

Dispensette S 4 page-2021

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BrandTech Scientific. Status as of: 0218/13

**For dispensing HF, we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring (Cat. No. 4640041).**



Toll Free: 888-522-2726  
www.brandtech.com