




MosID

identify. count. digitize.



Simple to use:



1. Load tray with mosquitoes

2. Insert tray into device

3. Capture image through website

4. Results available in seconds



Requirements:

- Flat 6in X 9in surface
- Wifi/Ethernet
- Computer
- AC wall power



Metrics:

- 98.3±0.7% species ID accuracy ¹
- 200-600 specimens per hour ²
- Streamlined data entry

¹ Preliminary data on 17 species, wild caught.

² Variable depending on whether pathogen testing is required.

Validation with partners through the 2021 season

Interested in being a Beta Customer? Take this one minute survey to let us know! Complete the survey to be entered into raffle for a \$20 Amazon gift card!

Survey link <https://tinyurl.com/vectech2021>