




# MosID

identify. count. digitize.



Simple to use:



1. Load tray with mosquitoes

2. Insert tray into device

3. Capture image through website

4. Results available in seconds



## Requirements:

- Flat 6in X 9in surface
- Wifi/Ethernet
- Computer
- AC wall power



## Metrics:

- 98.3±0.7% species ID accuracy<sup>1</sup>
- 200–600 specimens per hour<sup>2</sup>
- Streamlined data entry

<sup>1</sup> Preliminary data on 17 species, wild caught.

<sup>2</sup> Variable depending on whether pathogen testing is required.

## Validation with partners through the 2021 season

Interested in being a Beta Customer? Take this one minute survey to let us know! Complete the survey to be entered into raffle for a \$20 Amazon gift card!

Survey link <https://tinyurl.com/vectech2021>