



STERILE INJECTABLES • 503B PRODUCTS

503B VASOPRESSIN INJECTION BAGS FEE-FOR-SERVICE PROGRAM

You provide the vials, we prepare the bags.

WHAT IS IT?

It is a 503B Outsourced Compounding program — with a twist.

Our program is designed to offer customers the best of both worlds: meaningful savings purchasing vials at their lowest contracted price, paired with a guaranteed supply of premix, ready-to-use products made in accordance with cGMP.

It's producing dedicated batches specifically for dedicated customers.

HOW IT WORKS:

- The customer buys the Vasopressin[®] 10 mL vials through their usual channels who ship them directly to Athenex
- Athenex sends back to customers premix, ready-to-use 503B compounded 20 units/100 mL Vasopressin bags; made sterile-to-sterile
- For every one 10 mL vial Athenex receives, customers are guaranteed 10 bags in return
- Minimum order volume: 480 bags (*one full batch*)
- Customers determine the order volume and the order placement interval (*four weeks, six weeks or eight weeks*)
- Athenex and customer will develop schedule to optimize order placement, production and timely supply

BENEFITS:

Athenex Vasopressin Fee-For-Service Model

- Ready-to-use bags; guaranteed, just-in-time supply
- Manufactured in accordance with cGMP
- Premix, ready-to-administer products are recommended by ASHP*
- 90 days dating; room temperature; minimum waste
- Risk/cost of failed batch on Athenex
- Easy placement into automated dispensing stations

vs.


Hospital In-House (DIY)

- Vials are not ready to use; need to be admixed on site
- Potential to be admixed under less than cGMP conditions, such as USP <797>
- Time-sensitive drug; on-site admix brings potential for medication errors
- High potential for wastage
- Risk/cost of failed admix on hospital

*Whenever possible, medications should be available for inpatient use in unit-of-use and ready-to-administer packaging without further manipulation by the person administering the medication. ASHP Guidelines on Preventing Medication Errors in Hospitals; November, 2017¹

LIMITED RESOURCES?

Don't make it. Stock it!
Save time and labor!



Dating: 90 days | Room temp.

*****Sterile-to-Sterile Preparation*****

Athenex's premix, ready-to-use bags are made sterile-to-sterile and offer 90 days dating at room temperature.

